



**INNER POWER.
BETTER LIFE.**

SUMMARY


- Start the application p 3
- Create a client account p 4
- Body evaluations and treatment programs p 9
- Face evaluations and treatment programs p 19
- Photo assessment p 35
- Compare photo assessments p 40
- Send treatment program to your client p 43
- Send treatment program to your CELLU M6 ALLIANCE® p 46
- Free sessions p 51
- Backup p 53
- Updates p 56
- Data import (LPG SCAN tablet to the new one) p 60
- Data export p 62



**INNER POWER.
BETTER LIFE.**

START THE APPLICATION

Select the **endermoscan™** icon on your tablet's home screen to launch the application.






CREATE CLIENT ACCOUNT



INNER POWER.
BETTER LIFE.

CREATE CLIENT ACCOUNT

Once the application is open, select **NEW CLIENT** on the home screen.



INNER POWER.
BETTER LIFE.

CREATE CLIENT ACCOUNT

After discussing all contraindications with your client, please tick the **“Contraindications for endermologie® treatments were discussed during your evaluation”** box to start filling in the client’s profile.

LPG INNER POWER. BETTER LIFE. 14:38

HOME > CREATE AN ACCOUNT > CONTRAINDICATIONS

INFORMATION ON ENDERMOLOGIE® TREATMENTS

Strict general contraindications:

- Progressing cancer (except with written medical recommendation)
- Anticoagulant treatment
- Blood disease
- Disease with inflammatory eruption
- Infection or skin rash

Body endermologie® contraindications:

- Pregnancy (except for lower limbs)
- Phlebitis
- Organ transplant
- Kidney failure

Face endermologie® contraindications:

- Vitiligo
- Infectious, inflammatory acne
- Herpes

Localized Contraindications:

- Rosacea
- Varicosity and varicose veins
- Vein surgery (for 2 months)
- Areas that underwent plastic surgery (for 6 months)
- Skin lesions and wounds
- Keloid, hypertrophic, or recent scars
- Inguinal or abdominal hernias
- Lipoma
- Raised moles
- Filler injections (for 2 weeks)
- Botox® injections (for 2 weeks)
- Piercing
- Angioma
- Subcutaneous implants
- Gold and serrated thread lifts
- Acute osteo-articular or muscle injury

If in doubt, ask for a written medical opinion.

Contraindications for endermologie® treatments were discussed during your evaluation.

CREATE CLIENT ACCOUNT

- Fill in client information.
- Click on **CREATE THIS CLIENT PROFILE.**
- THE BMI is calculated automatically after filling in the height and the weight.

The screenshot shows a mobile application interface for creating a client account. At the top, the LPG logo and tagline 'INNER POWER. BETTER LIFE.' are visible, along with the time 14:39 and a Wi-Fi signal icon. Below the header, a breadcrumb trail reads 'HOME > CREATE AN ACCOUNT > PERSONAL INFORMATION'. The main form area contains the following elements:



- Gender selection: Woman Man
- Text input fields for 'Last Name *' and 'First Name *'.
- A text input field for 'Date of Birth *'.
- Text input fields for 'Address', 'City', and 'Zip Code'.
- Text input fields for 'Phone' and 'Email'.
- Text input fields for 'Height (in cm)' and 'Weight (in kg)'.
- Text input fields for 'medical follow-up' and 'surgical procedures'.
- A profile picture placeholder with a pencil icon for editing.
- A back arrow on the left side.
- A large button at the bottom labeled 'CREATE THIS CLIENT PROFILE'.



INNER POWER.
BETTER LIFE.

CREATE CLIENT ACCOUNT

- To view or to edit a client profile, select your client's name and click on 
- To search for a client in the existing database, type the client's name in the **SEARCH** bar.
- To archive a client profile, click on 
- To delete a client profile, click on 






BODY EVALUATION & TREATMENT PROGRAMS

LPG

INNER POWER.
BETTER LIFE.

BODY EVALUATION

Once the client account has been created, choose to carry out a **body evaluation** by selecting the concerned area.



BODY EVALUATION

It is important to inform your client of the different endermologie® treatment programs success levers.

The screenshot shows a mobile application interface. At the top, there is a dark grey header bar with the LPG logo on the left, the text "INNER POWER. BETTER LIFE." in the center, and the time "14:42" on the right. Below the header is a white navigation bar with a hamburger menu icon on the left, the text "HOME > LAURE DURAND > NEW" in the center, and a Wi-Fi signal icon on the right. The main content area is a light grey slide with a white background. The slide is titled "THE WINNING TRIO" in bold black text. Below the title, there is a paragraph: "Results will depend on her/his ability to combine endermologie® treatments with a balanced diet and regular exercise." This is followed by two sections: "A BALANCED DIET" and "REGULAR EXERCISE", each with a list of bullet points. The slide has a left arrow and a right arrow on its sides, indicating it is part of a carousel.

LPG INNER POWER. BETTER LIFE. 14:42

HOME > LAURE DURAND > NEW

THE WINNING TRIO

Results will depend on her/his ability to combine endermologie® treatments with a balanced diet and regular exercise.

A BALANCED DIET

- 5 fruits and vegetables per day; proteins, carbohydrates, and various foods at each meal; limit sweet, salty, and fatty foods.
- It is recommended that you drink regularly to keep your body hydrated (ideally between meals).

REGULAR EXERCISE

- The equivalent of walking at least 30 minutes per day is recommended.
- For fat-elimination results, it is recommended that you exercise after your body endermologie® session.

BODY EVALUATION

It is important to inform your client of the different endermologie® treatment programs success levels. (i.e.: Results may vary depending on to the hormonal stability)


The screenshot shows a mobile application interface for a 'FEMALE EVENTS' evaluation. At the top, there is a dark grey header with the 'LPG' logo, the tagline 'INNER POWER. BETTER LIFE.', the time '14:44', and a Wi-Fi signal icon. Below the header is a navigation bar with a hamburger menu icon, the text 'HOME > LAURE DURAND > NEW', and a right-pointing arrow. The main content area is titled 'FEMALE EVENTS' and contains three sections: 'PREGNANCY' with radio buttons for 'none', '1', '2', '3', and '3+' (where 'none' is selected); 'CONTRACEPTION' with radio buttons for 'no' and 'yes' (where 'yes' is selected); and 'MENOPAUSE' with radio buttons for 'no' and 'yes' (where 'no' is selected). Below these sections is a 'COMMENTS' field. Navigation arrows are visible on the left and right sides of the form area.



INNER POWER.
BETTER LIFE.

BODY EVALUATION

- Fill in client's desired objectives.
- Select the type of treatment desired according to your client's objectives :
 - Slimming (*endermologie® ID treatments*)
 - Well-being treatments (*these do not require evaluation*)




The screenshot shows the LPG mobile application interface. At the top, the LPG logo and tagline "INNER POWER. BETTER LIFE." are on the left, the time "14:44" is in the center, and a Wi-Fi icon is on the right. Below the header, a navigation bar shows "HOME > LAURE DURAND > NEW". The main content area features a text input field labeled "Objective:" with the placeholder text "Fill in desired objectives here". Below this field are two selectable options: "Slimming", which is accompanied by an image of a woman's midsection and a hand icon pointing to it, and "Well-being", which is accompanied by an image of a woman performing a backbend yoga pose. A blue hand icon is shown clicking on the "Slimming" option.

BODY EVALUATION

SLIMMING


- Take 4 photos (front, right profile, left profile, back). Maintain you tablet in vertical (portrait) format. The photo must be taken when level sensors are green. Report to the integrated markers of the photo assistant to position your client properly. Once the picture is taken, click on **REGISTER** or **CANCEL** if you are not satisfied.
- To take other pictures, click on **ADD ADDITIONAL PICTURE** and choose a title.



BODY EVALUATION

Water retention questionnaire

Fill it in if your client is concerned with this issue to determine the endermologie® cosmetics best adapted to their needs.



The screenshot shows a mobile application interface for a 'Water retention questionnaire'. The app's header includes the LPG logo and the tagline 'INNER POWER. BETTER LIFE.' on the left, the time '14:44' in the center, and a Wi-Fi signal icon on the right. Below the header is a navigation bar with a hamburger menu icon, the text 'HOME > LAURE DURAND > NEW', and a right-pointing arrow. The main content area is titled 'Water retention questionnaire' and contains five questions, each with multiple-choice options:


- 1 - Where is your cellulite?
 - diffuse
 - in one area
- 2 - Can your weight vary 1 kg to 3 kg within a 48-hour period?
 - yes
 - no
- 3 - How do your legs feel?
 - Only heavy and painful at the end of the day
 - Already swollen in the morning and worsening at the end of the day
 - Only swollen during my premenstrual period
 - No problem with my legs
- 4 - Describe your ankles.
 - They are swollen after prolonged sitting (for example, a flight)
 - They are always swollen with edema around the ankle bones
 - Nothing in particular regarding my ankles
- 5 - Do you have any right now?
 - Swollen hands and feet?
 - Tree-trunk legs?
 - Noticed a varying increase in volume from one week to the next?
 - A cortisone-based medical treatment?



INNER POWER.
BETTER LIFE.

BODY EVALUATION

- Select the client's targeted areas.
- Select the type of tissue (supple, dense, fibrous) corresponding to each selected area.
- Adjust treatment time corresponding to each area.



BODY EVALUATION

Recommended endermologie® protocols

Endermologie® ID treatment(s) previously selected are automatically checked.
According to your client's objectives, you can add a well-being treatment.

The screenshot displays the LPG mobile application interface. At the top, the LPG logo and the slogan "INNER POWER. BETTER LIFE." are visible on the left, and the time "14:45" and a Wi-Fi icon are on the right. Below the header, a navigation bar shows "HOME > LAURE DURAND > NEW". The main content area is titled "Recommended endermologie® protocols" and is divided into three sections:



- endermologie® ID treatment(s)**
 - Back
 - Stomach
 - Back of thighs
 - Inner thighs
- Specific treatment(s)**
 - Young mum
- Well-being treatment(s)**
 - Relaxation
 - Fluid Mobilization
 - Endermopuncture

Navigation arrows are present on the left and right sides of the main content area.

BODY EVALUATION

Final evaluation and protocol, cosmetics & endermonutri recommendation

- Adjust the number of recommended sessions.
- To confirm the evaluation and to create the customized program, select **SUBMIT**.
- Your client's program is then saved.





FACE EVALUATION & TREATMENT PROGRAMS

LPG

**INNER POWER.
BETTER LIFE.**

FACE EVALUATION

To start a new evaluation, click on



LPG INNER POWER. BETTER LIFE. 14:45

HOME > LAURE DURAND

BODY 14 JUNE 2024

- Back - Stomach - Back of thighs - Inner thighs
- Relaxation

Two photos of a person in white underwear, a grid icon, a camera icon, and a collection of LPG skincare products.

Share, Edit, Delete icons.


Bottom navigation bar with a plus sign (+) and a hand cursor pointing to it.



INNER POWER.
BETTER LIFE.

FACE EVALUATION

To start a face evaluation, click on the face area.



INNER POWER.
BETTER LIFE.

FACE EVALUATION

Fill in client's desired objectives.


The screenshot shows a mobile application interface. At the top, there is a dark grey header bar with the 'LPG' logo and the tagline 'INNER POWER. BETTER LIFE.' on the left, the time '14:46' in the center, and a Wi-Fi signal icon on the right. Below the header is a white navigation bar with a hamburger menu icon on the left and the breadcrumb 'HOME > LAURE DURAND > NEW' in the center. The main content area is a large white rectangle with a thin border. Inside this area, the text 'Objective:' is at the top left, followed by a larger text input field containing the placeholder text 'Fill in desired objectives here'. On the left and right sides of the main content area, there are grey vertical bars with black arrow icons pointing left and right, respectively, indicating a scrollable view.



INNER POWER.
BETTER LIFE.

FACE EVALUATION

- Take 3 photos (right profile, face, left profile). Maintain you tablet in vertical (portrait) format. The photo must be taken when level sensors are green. Report to the integrated markers of the photo assistant to position your client properly. Once the picture is taken, click on **REGISTER** or **CANCEL** if you are not satisfied.
- To take other pictures, click on **ADD ADDITIONAL PICTURE** and choose a title.



INNER POWER.
BETTER LIFE.

FACE EVALUATION

Select the corresponding skin aging type.

LPG INNER POWER. BETTER LIFE. 14:46

HOME > LAURE DURAND > NEW

Choose the type of aging

< >

Sagging Hollowing Thickening



INNER POWER.
BETTER LIFE.


FACE EVALUATION

Select up to 4 priority areas.
For each selected area, choose the main objective.

The screenshot displays the LPG mobile application's face evaluation screen. At the top, the LPG logo and tagline 'INNER POWER. BETTER LIFE.' are visible, along with the time '14:46' and a Wi-Fi icon. Below this is a navigation bar with a menu icon and the text 'HOME > LAURE DURAND > NEW'. The main content area features a portrait of a woman with several teal selection boxes overlaid on her face and neck. Lines connect these boxes to labels: 'Wrinkles and Fine Lines 10 min' points to the eye area, 'Blotchy' points to the right cheek, 'Stored Fat 10 min' points to the chin, and 'Wrinkles 10 min' points to the neck. A large '< 30 min' is positioned on the left side of the image, and a right-pointing arrow is on the right side. The bottom of the screen shows a horizontal line and a vertical line forming a corner.

FACE EVALUATION

Move the cursor to select the corresponding level.



INNER POWER.
BETTER LIFE.

FACE EVALUATION

Move the cursor to select the corresponding level.

The screenshot shows a mobile application interface for a face evaluation tool. At the top, there is a dark grey header with the 'LPG' logo and the slogan 'INNER POWER. BETTER LIFE.' on the left, the time '14:47' in the center, and a Wi-Fi signal icon on the right. Below the header is a white navigation bar with a hamburger menu icon on the left and the breadcrumb 'HOME > LAURE DURAND > NEW' in the center. The main content area is titled 'CHIN' and features a profile photograph of a person's chin. Below the photo is a horizontal scale with six numbered points (1 to 6). A teal line highlights the segment between points 1 and 3. A hand icon is positioned over point 3, indicating it is the selected level. Navigation arrows are visible on either side of the scale.

FACE EVALUATION

Move the cursor to select the corresponding level.

The screenshot displays the LPG mobile application interface. At the top, the LPG logo and the slogan "INNER POWER. BETTER LIFE." are visible on the left, the time "14:47" is in the center, and a Wi-Fi icon is on the right. Below this is a navigation bar with a menu icon, the text "HOME > LAURE DURAND > NEW", and a right-pointing arrow. The main content area is titled "NECK" and features a photograph of a person's neck. Below the photo is a horizontal scale with six numbered points (1 to 6). A hand cursor is positioned over the number 3, indicating the selected evaluation level. The scale is highlighted with a teal line from point 1 to point 3.



INNER POWER.
BETTER LIFE.

FACE EVALUATION

Select most adapted protocol(s) according to your client's objectives.

The screenshot displays the LPG mobile application interface. At the top, the LPG logo and the slogan "INNER POWER. BETTER LIFE." are on the left, and the time "14:50" and a Wi-Fi icon are on the right. Below the header, a navigation bar shows "HOME > LAURE DURAND > NEW". The main content area is titled "Recommended endermologie® protocols" and is divided into three sections: "Targeted treatment(s)", "Classic treatment(s)", and "Signature treatment(s)".

Recommended endermologie® protocols

Targeted treatment(s)

- Eye Wrinkles
- Double Chin
- Neck

Classic treatment(s)

- Fluid Mobilization
- Age-defying Replumping
- Age-defying Resculpting
- Age-defying Skin Toning

Signature treatment(s)

- Perfect Eyes & Lips





INNER POWER.
BETTER LIFE.

FACE EVALUATION

Final evaluation and protocol, cosmetics & endermonutri recommendation

- Adjust the number of recommended sessions.
- To confirm the evaluation and program creation, select **SUBMIT**.
- Your client's program is then saved.



INNER POWER.
BETTER LIFE.

VIEWING PROGRAMS

SESSIONS

To view the completed endermologie® sessions (protocol description, date). This will enable a precise treatment program follow-up.

PHOTOS ASSESSMENT

To make, to consult and to compare photo assessments.

COSMETICS


To view cosmetics associated with your client's programs.

ENDERMONUTRI

To view food supplements associated with your client's programs.

List of programs

Date of the first evaluation



To create new facial or body treatment program

To send your client's program by email

To modify your client's program

To delete your client's program



INNER POWER.
BETTER LIFE.

VIEWING SESSIONS

 : To confirm a session in case of connection loss with your CELLU M6 ALLIANCE®

 : To comment a session (settings used...)

 : To delete a session

 : To add a session



The screenshot shows the LPG mobile app interface. At the top, the status bar displays 'LPG INNER POWER. BETTER LIFE.', the time '14:50', and a Wi-Fi signal icon. Below the status bar, a navigation bar shows 'HOME > LAURE DURAND > MY SESSIONS'. The main content area displays a list of sessions under the heading 'BACK - STOMACH - BACK OF THIGHS - INNER THIGHS' with a progress indicator '0/10 40 min' and an add button '+'. Each session entry includes a session number (1/10, 2/10, 3/10), a list of exercises with their durations, and three action icons: a calendar icon, a comment icon, and a delete icon.

Session	Exercise	Duration	Actions
1/10	Back - dense	10 min	Calendar, Comment, Delete
	Stomach - supple	10 min	
	Back of thighs - dense	15 min	
	Inner thighs - dense	5 min	
2/10	Back - dense	10 min	Calendar, Comment, Delete
	Stomach - supple	10 min	
	Back of thighs - dense	15 min	
	Inner thighs - dense	5 min	
3/10	Back - dense	10 min	Calendar, Comment, Delete
	Stomach - supple	10 min	
	Back of thighs - dense	15 min	
	Inner thighs - dense	5 min	
	Back - dense	10 min	

IN EXPERT MODE

*Possibility to save the settings used during the session for next sessions. *Possibility to save the time spent on each step.



INNER POWER.
BETTER LIFE.

VIEWING COSMETICS

You can view cosmetics associated with each treatment program.

The screenshot shows the LPG mobile application interface. At the top, the header includes the LPG logo with the tagline "INNER POWER. BETTER LIFE.", the time "10:47", and a Wi-Fi signal icon. Below the header, a navigation bar shows "HOME > LAURE DURAND > MY BEAUTY ROUTINE".

The main content area is split into two columns. The left column is a menu with several items, each with a small product image and a title:

- Eye Wrinkles - Double Chin - Neck
- Perfect Eyes & Lips
- Preparing Micellar Water** (highlighted in a teal bar)
- Post-Care Eye Contour Mask
- Eye & Lip Contour Cream
- Anti-Aging Renewal Serum
- Anti-Aging Renewal Cream

The right column displays the details for the selected product, "PREPARING MICELLAR WATER". It features a larger product image, a list of benefits, and usage instructions:

- Removes impurities and makeup Moisturizes and soothes

Soak a cotton pad and cleanse face, eyes and neck, morning and evening. Dry the skin gently without rinsing.

VIEWING ENDERMONUTRI

You can view food supplements associated with each treatment program.

The screenshot shows the LPG mobile app interface. At the top, the LPG logo and tagline 'INNER POWER. BETTER LIFE.' are on the left, the time '15:37' is in the center, and a Wi-Fi icon is on the right. Below this is a navigation bar with a hamburger menu icon, 'HOME > LAURE DURAND > MY BEAUTY ROUTINE'. The main content area is divided into two columns. The left column is a list of supplements under the heading 'EYE WRINKLES - DOUBLE CHIN - NECK'. The items are: 'VITALITY BOOSTER' (highlighted with a teal background), 'HYALURONIC ACID', 'COLLAGEN', 'OMEGA', and 'WATERGLOW'. The right column is titled 'VITALITY BOOSTER' and features a product image of a glass jar with a metal lid. To the right of the image is a bulleted list of benefits: 'Stimulates vitality and physical & mental performance', 'Reduces fatigue', and 'Strengthens the immune system'. Below the image and list is a paragraph of instructions: '2 capsules per day preferably in the morning and before 3 pm to prevent any difficulty in falling asleep. For optimal results, it is recommended to follow this program for at least 2 consecutive months. It may be repeated at the start of any endermologie® program and/or change in season.'




PHOTO ASSESSMENT

LPG

INNER POWER.
BETTER LIFE.

FACE PHOTO ASSESSMENT

After clicking on the “Photo assessment” icon, click on  to make a new photo assessment. All photos taken during the first evaluation are saved here.



 : To save photos assessment on your tablet

 : To delete photos assessment



INNER POWER.
BETTER LIFE.


FACE PHOTO ASSESSMENT

The first evaluation's objectives are mentioned (reminder).
You can add comments (number of sessions, results..)

The screenshot shows a mobile application interface. At the top, there is a dark grey header bar with the 'LPG' logo on the left, the text 'INNER POWER. BETTER LIFE.' in the center, the time '15:55' on the right, and a Wi-Fi signal icon. Below the header is a white navigation bar with a hamburger menu icon on the left and the breadcrumb path 'HOME > LAURE DURAND > PHOTO ASSESSMENT'. The main content area is a light grey rectangle containing two text input fields. The first field is labeled 'Objective:' and is empty. The second field is larger and contains the placeholder text 'Add your comments here'. To the left and right of this second field are grey arrow icons pointing left and right, respectively, indicating a scrollable area.


FACE PHOTO ASSESSMENT

- Take 3 photos (right profile, face, left profile). Maintain you tablet in vertical (portrait) format. The photo must be taken when level sensors are green. Report to the integrated markers of the photo assistant to position your client properly. Once the picture is taken, click on **REGISTER** or **CANCEL** if you are not satisfied.
- To take other pictures, click on **ADD ADDITIONAL PICTURE** and choose a title.



FACE PHOTO ASSESSMENT

Move the cursor to select the corresponding level.



INNER POWER.
BETTER LIFE.




COMPARE PHOTO ASSESSMENTS



INNER POWER.
BETTER LIFE.

COMPARE FACE PHOTO ASSESSMENTS


Click on 



The screenshot shows the LPG mobile app interface. At the top, the status bar displays the time 14:54 and a Wi-Fi signal icon. Below the status bar, the app header includes the LPG logo and the slogan "INNER POWER. BETTER LIFE." followed by a navigation breadcrumb: "HOME > LAURE DURAND > PHOTO ASSESSMENT". The main content area features two identical rows, each representing a photo assessment. Each row contains a profile photo of a woman, a magnifying glass icon over a document, and the date "14 JUNE 2024". To the right of each row, it says "0 PHOTOS" and includes a right-pointing chevron, a download icon, and a trash icon. At the bottom of the screen, there is a dark grey navigation bar with three icons: a plus sign, a compare icon, and a hand icon pointing to the compare icon.

COMPARE FACE PHOTO ASSESSMENTS

- Using the arrows, select the photo assessments
- Next, select the picture: face, back, left profile, right, profile
- To save the comparison, click on 
- Choose a title and save it in a folder on your tablet.






SEND TREATMENT PROGRAM TO YOUR CLIENT



INNER POWER.
BETTER LIFE.

SEND TREATMENT PROGRAMS


- Choose the concerned program and click on 
- Select cosmetics and food supplements to import in the document. Next, click on **OK**.



INNER POWER.
BETTER LIFE.

SEND TREATMENT PROGRAM

After setting up your email account, customize the message which will be sent.
The program (PDF document) is automatically enclosed with your email. Next, click on **SEND**.



INNER POWER.
BETTER LIFE.



SEND TREATMENT PROGRAM
TO YOUR CELLU M6 ALLIANCE®


LPG

INNER POWER.
BETTER LIFE.

TREATMENT PROGRAM LAUNCH

After connecting your tablet to your CELLU M6® device (please refer to the Tablet Quickstart guide if necessary), select **CUSTOMIZED SESSION**.


If your device is connected to your tablet, this icon  appears on your tablet's screen.



INNER POWER.
BETTER LIFE.

TREATMENT PROGRAM LAUNCH

Select the concerned client in your database.



Your client's programs can be downloaded when the tablet is connected. They will be accessible at any time. The results will be sent as soon as the tablet reappears.



**INNER POWER.
BETTER LIFE.**

TREATMENT PROGRAM LAUNCH

Select the desired program.



In case of connection loss, possibility to update manually the sessions on the machine validated on the tablet by clicking on 



INNER POWER.
BETTER LIFE.

TREATMENT PROGRAM LAUNCH

The session can start!

LPG INNER POWER. BETTER LIFE. 15:34

LAURE DURANT > STOMACH - BUTTOCKS - SADDLEBAGS > STOMACH > DENSE 00:00

Step 1/22 4:35 -5:25

RESURFACING CREAM

Before putting on LPG® endermowear™, apply the cream on the entire body by massaging upwards until completely absorbed giving special attention to stretch marks areas.



The image shows three white bottles of LPG skincare products. The two bottles in the foreground are 'CRÈME MICRO-PEELING RESURFACING CREAM' with a volume of 400 ML e / 13.5 FL. OZ. The bottle in the background is 'PRÉPARATEUR DE SOIN CORPS SKINCARE BODY PREPARER' with a volume of 400 ML e / 13.5 FL. OZ. All bottles feature the LPG logo and list ingredients: AHA, Antioxydants, Vitamines.



INNER POWER.
BETTER LIFE.




FREE SESSIONS



INNER POWER.
BETTER LIFE.

VIEWING FREE SESSIONS

To access your client's free sessions, click on 
Next, click on **MY FREE SESSIONS**.



The screenshot shows the LPG mobile app interface. At the top, the header includes the LPG logo, the slogan "INNER POWER. BETTER LIFE.", the time "16:12", and a Wi-Fi signal icon. Below the header, a navigation bar shows "HOME > LAURE DURAND". The main content area displays two session cards for "FACE" and "BODY", both dated "14 JUNE 2024".

FACE (14 JUNE 2024):

- Eye Wrinkles - Double Chin - Neck
- Perfect Eyes & Lips

BODY (14 JUNE 2024):

- Back - Stomach - Back of thighs - Inner thighs
- Relaxation

Each session card includes a front and back view of the client, a calendar icon, a camera icon, a product image, and a set of action icons (share, edit, delete). A blue hand icon points to a menu icon in the top right corner of the session cards, which has opened a dropdown menu with the following options:

- My free sessions
- Imported evaluations



BACKUP



INNER POWER.
BETTER LIFE.

BACKUP





Should you skip this essential step and should any problem occur (operating error, technical failure, theft...), you will lose your data

For more security, save all your data on an external support of your choice: micro SD card available in supermarkets or in specialized stores, personal computer, external hard drive...

By default, your database will be stored on your tablet.

Click on the menu button on the upper left side of the screen and select **SETTINGS**. Click on **DATABASE**.

Click on « **BACKUP DATABASE** » and choose your micro SD card as backup location. The number of possible backups is conditioned by your micro SD card's memory size.



**INNER POWER.
BETTER LIFE.**

AUTOMATIC BACKUP

To plan an automatic backup, click on **AUTOMATIC BACKUP**.



Click on the cursor to activate the automatic backup



Select the periodicity, the location.

Enter your password for security protection. Please keep it safe because without it, it will be impossible to restore your backups. For privacy reasons, LPG SYSTEMS does not have access to your password and cannot ask you for it. LPG SYSTEMS can not be held responsible for any loss of your data.

Click on **OK**.



INNER POWER.
BETTER LIFE.




UPDATES



INNER POWER.
BETTER LIFE.



UPDATES

When an update is available, the following screen appears when you start the machine.
The "i" allows to get information about the update.
Click on **OK** to start the update.



UPDATES

Please follow the instructions on the screen and continue by clicking on **OK**.



Enter your password for security protection. Please keep it safe because without it, it will be impossible to restore your backups. For privacy reasons, LPG SYSTEMS does not have access to your password and cannot ask you for it.



LPG SYSTEMS can not be held responsible for any loss of your data.



INNER POWER.
BETTER LIFE.

UPDATES

At the end of the backup, your tablet can display this message. Click on **Settings**. Allow the installation of unknown sources with the cursor.



INNER POWER.
BETTER LIFE.



DATA IMPORT (LPG® SCAN tablet to the new one)

LPG



INNER POWER.
BETTER LIFE.

LPG® SCAN DATABASE

Before importing, please save your database.

To import your LPG SCAN database (client profile, photos), please insert a micro SD card into the intended location on the tablet. Then click on the menu button on the upper left side of the screen and select **SETTINGS**.

In **DATABASE**, click on **IMPORT DATABASE** and choose a backup location.





DATA EXPORT





INNER POWER.
BETTER LIFE.

DATA EXPORT

To access your tablet's settings, click on 

« **Disclaimer** » makes you aware of confidentiality and integrity of data if you decide to export our customer data (please check *the export of customer data*).



LPG® is not responsible for the use of data outside of its application.