

Registered trademark

OPERATOR'S MANUAL HUBER 360® FIIT

Please read this entire manual carefully before using your device.

Keep this manual and the separate installation instructions for future reference

 $\hfill LPG$ Systems 2024. LPG and HUBER are registered trademarks. Reproduction, even in part, is strictly prohibited.

GU1401 - GB Indice C1 of 01/24



Congratulations on the purchase of your HUBER 360[®] FIIT device.

This model represents the culmination of many years of research, combining the latest technologies with LPG's expertise.

You can fully appreciate the technical perfection and reliability that have made LPG Systems the leader in this field.

This operator's manual describes how to operate the device, instructions for periodic maintenance, and safety guidelines.

With the aim to continuously provide satisfaction, your device is equipped with a specific software devised to ensure connection with the LPG dedicated server. The data collected through this software will allow LPG to improve its support and maintenance servcie.

Your device should be used only by a professional who has completed the manufacturer's training provided by LPG Systems or an approved provider, if you live outside of France.

If you have any questions about operating or maintaining your device, please contact LPG via the Customer Service Department of your local LPG Distributor.

7 WARNING

To better respond to customer requirements and expectations, LPG Systems is constantly looking for ways to improve the design and quality of its products. For this reason, there may be minor differences between your equipment and the unit described in this manual.

\rightarrow **PACKAGE CONTENTS**

4


- → One HUBER 360[®] FIIT
- \rightarrow One heart rate monitor
- \rightarrow One tablet
- \rightarrow One operator's manual
- \rightarrow One power cable
- \rightarrow One set of quick unpacking and set-up instructions
- \rightarrow One POS display



\rightarrow contents


<u>1</u> .	DEVICE DESCRIPTION	7
<u>2</u> .	CONTROLS AND ACCESSORIES	9
	 2.1 Control Screen. 2.2 Handles. 2.6 Heart rate monitor. 2.7 Tablet 	
3.	IMPORTANT SAFETY INFORMATION	
	 3.1 Warning. 3.2 Precautions. 3.3 Clearance. 3.4 Warnings. 	
<u>4</u> .	MAINTENANCE	
	4.1 Cleaning the device.4.2 Identification nameplate.4.3 Power cord.	
5.	TROUBLESHOOTING	
6.	TECHNICAL SPECIFICATIONS	
<u>7</u> .	WARRANTY	
	TERMS AND CONDITIONS OF THE WARRANTY	
	LIMITATIONS AND EXCLUSIONS	

\rightarrow DEVICE DESCRIPTION



↗ WARNING

The unit can only operate if its power cable is connected to the power outlet and if the on/off switch is turned on.




ACCESSORIES | CONTROLS AND

9

\rightarrow handles



\rightarrow **PLATFORM**



\rightarrow HEART RATE MONITOR

Your HUBER 360 $^{\otimes}$ FIIT may come with a Heart Rate Monitor that constantly displays your heart rate.

Activate your Heart Rate Monitor by centering it on your chest.

Moisten the contact areas with water or gel. Next, follow the instructions on the screen.

Your heart rate will appear at the lower right corner of the screen.



Heart Rate

- · Display: 30 to 240 bpm
- · Limit: 30 to 230 bpm
- · Maximum percentage for display: 240 bpm

✓ WARNING

The heart rate meter is designed only for casual use. It should not be used for obtaining measurements that require professional or laboratory-quality precision.

The serial number for the heart rate meter can be found on the back of the belt.

The heart rate monitor contains a battery that cannot be disposed of with household garbage, it requires special recycling. When the battery dies, dispose of it in an authorized location for recycling.

\rightarrow TABLET

Your HUBER 360 $^{\circ}$ FIIT may come with a tablet with which you can download and access the HUBER 360 $^{\circ}$ FIIT application.

The HUBER 360 $^{\circ}$ FIIT application is downloadable thanks to the "my LPG" application, available on the Play Store.

This application enables you to:

- carry out an assessment
- create customized programmes
- track progression
- monitor client at a distance
- send programmes to clients

Details about how to use the application are provided when you open it.



12

→ IMPORTANT SAFETY INFORMATION

The following safety precautions must always be observed when using an electrical device. Read all instructions before using the device.

DANGER - TO MINIMIZE THE RISK OF ELECTRIC SHOCK:

- Always unplug the device from the power outlet after use and before cleaning it (ensure that the power switch is off).
- Check that the voltage indicated on the nameplate matches the voltage of the power outlet.
- The unit must be connected to a grounded outlet using the supplied power Cord¹, in accordance with current electrical standards. Electrical adapters must not be used with this equipment.
- Ensure that a differencial device (SI type, detecting the DC and AC components from a frequency inverter) is installed on the circuit upstream from the machine. This must comply strictly with all appropriate safety regulations and restrictions.
- Ensure that the circuit-breaker rating is in correct correspondence to the HUBER 360[®] FIIT (the upstream circuit breaker must be rated higher than the machine).
- The unit must be connected to a dedicated single power line.

\rightarrow WARNING

TO MINIMIZE THE RISK OF BURNS, FIRE, ELECTRIC SHOCK, OR PERSONAL INJURY:

- The device should never be left unattended while plugged in.
- Unplug the device when not in use for a long period of time or before installing or removing parts.
- Carefully monitor the device when it is used by, on, or around children or disabled individuals.
- Children should be supervised to ensure that they do not play with the appliance.
- This device is not intended for use by persons (including children) whose physical, sensory or mental capabilities are reduced, or persons without experience or knowledge, unless they have benefited from the intermediary of a person responsible for their safety, surveillance or prior instructions concerning the use of the device.
- Never use this device if the power cord¹ or outlet is damaged, if it is not working properly, if it has been dropped, or if it has been exposed to excessive moisture.
- Return the device to the LPG Systems service center for examination and repair.
- Never use the device if the vents are blocked. Keep openings free of dust and other contaminants.
- Do not drop or insert any objects into any of the openings. This could damage the machine.
- Do not use the device on a dirt or uneven floor or in a humid or salty environment, and do not expose it to bad weather.
- Cleaning and maintenance of the device by the user should not be performed by unattended children.
- Do not use the device around aerosol (sprays) or oxygen products.

14

¹ Provided power cord information: Europe VII-H05VVF3G 2,50-g25

▼ WARNING

→ PRECAUTIONS

WARNING : KEEP THESE INSTRUCTIONS

- As with any physical training, a medical examination should be performed before using the equipment.
- The professional should ensure that the settings (platform oscillation speed, amplitude, activity/rest duration, etc.) are correct for the user's morphology and physical condition.
- The equipment should be used under careful supervision by a trained professional.
- Users should warm up properly before using the equipment for exercise.
- Users should breathe regularly while exercising.
- Using this equipment with maximum effort is not recommended while the platform is moving.
- The user should stop immediately if experiencing any discomfort during exercise.
- Do not operate the unit in unsuitable environmental conditions (see technical specifications).
- Do not load the machine with more than the recommended weight (see technical specifications)
- Do not move the equipment without an LPG technician being present.
- As its name implies, the emergency shutoff button should be used only in emergencies.
- The tablet should only be used by a trained professional near the machine.
- Do not use the tablet and the user interface simultaneously when controlling the machine.
- This symbol means "High Voltage":

\rightarrow CLEARANCE

Provide a clearance of 1 meter all the way around the device, keeping the area free of any objects or people.



\rightarrow WARNINGS

- The equipment should not be used by anyone with heart, respiratory, neurological, or rheumatic disease who is advised against physical activity.
- In the case of regaining fitness, maintenance and preventative activity, the patient should be medically certified to undertake physical activity.
- Due to the risk of possible interference, it is important that the professional ensure that the patient is not fitted with a personal medical equipment, such as a pacemaker or a hearing aid.
- If this is the case, determine with certainty whether the equipment in question emits a signal that may interfere with the functioning of the unit.
- The unit must not be used by someone with a medical condition unless accompanied at all times by a qualified professional certified to the standards of the country of use.

Z WARNING


This machine contains programs meant to guide the user in order for them to gain the best possible results that can be expected in their case. This in no way guarantees that the care will be successful. Success varies according to each user's morphology and physiology. It is strictly prohibited to modify your unit without the authorization of LPG.

\rightarrow CLEANING THE DEVICE

You should clean your device as often as possible, not only for reasons of hygiene and attractiveness, but also because cleaning helps maintain it in good working order and prolong its life:

- Use a damp sponge to clean the outside of all covers.
- Use wipes impregnated with a bactericidal and fungicidal solution to clean the handles.
- Use a vacuum cleaner with a narrow attachment to clean the platform and the inside of the accessory positioning recesses.

→ IDENTIFICATION NAMEPLATE



18

▼ WARNING

The identification plate could be modified. The identification plate apposed under your device is the approved identification plate.

→ IDENTIFICATION NAMEPLATE (CONT.)

The identification nameplate is located on the back of your device.

The nameplate indicates the voltage of the device and the serial number to identify your device. When contacting LPG Systems customer service for a

technical problem, please provide the serial number for your HUBER 360[®] FIIT.

This serial number gives you the year and month of manufacture of your device.

The letter indicates the year of manufacture: A=2010, B=2011..., E=2014,...


The two digits indicate the production month: 01=January, 02=February,

This symbol means that the unit was sold after August 13, 2006. In conformity with the 2002/96/EC directive, it cannot be thrown away with standard household waste but must be disposed of in an appropriate sorted collection. When your device reaches the end of its useful life, it must be brought to an appropriate recycling center or returned to your dealer. By doing so, you help the environment by contributing to the conservation of natural resources and the protection of human health.

- This symbol means that some specific warnings or precautions associated with this device are not on the label.
- This symbol means, "Refer to the operator's manual".

This symbol means that the the maximum patient weight permitted on this equipment is 140Ka

This symbol means store the device somewhere protected from inclement weather.



- This symbol indicates relative humidity limits.
- This symbol means "Danger: High Voltage"

\rightarrow power cord

If the power cord of your device is damaged, please contact LPG Systems Customer Service via the Customer Service Department of your local LPG distributor, for a replacement.

NOTES

\rightarrow TROUBLESHOOTING

If your unit is not working properly, please run through the following checklist before calling LPG Systems Customer Service:

- Is the unit properly connected to a power outlet?
- Is the power outlet live?
- Is the power switch in the ON position?
- Is the emergency stop button pressed? If so, turn off the power, release the emergency stop button, and then restart the machine.
- For the symptoms listed below, stop the machine and wait 60 seconds before restarting:
- Error message displayed on the screen;
- Screen not working;
- Excessive weight applied to the platform during exercise;
- Fault on the sensor.

Once these checks have been completed and if the problem persists, please contact LPG Systems Customer Service or the nearest authorized dealer, indicating the model of your unit and its serial number.

\rightarrow TECHNICAL SPECIFICATIONS

• General Features

Dimensions l x w x h	cm
Net weight	kg
Maximum weight allowed on the platform	kg

• Electrical Features

Voltage	V
Rated frequency	z
Power consumption 1800VA - Momentary 1800VA	4

• Platform Movement

Maximum oscillation frequency 1	tr/sec
Maximum amplitude	10°

• Force Measurement

Handle measurement scale	0 to 85 kg
--------------------------	------------

• Weight measurement

Measurement range.....0 to 165 kg

User Replaceable Batteries

Heart rate monitor	3V CR2032 battery
--------------------	-------------------

• Surrounding environment restrictions

Operating temperature+ 15°C to + 40°	°C
Storage temperature 20°C to + 70°	°C

Ambient relative humidity	
Atmospheric pressure	
Maximum altitude	

$\rightarrow~$ TECHNICAL SPECIFICATIONS (CONT.)

Device designed for uninterrupted use.

HUBER 360[®] FIIT is **C** marked in accordance with applicable legislation and standards under the emc directive 2004/108/ec, low voltage directive 2006/95/EC and 2006/35/UE machine.

ELECTROMAGNETIC COMPATIBILITY

- HUBER 360[®] FIIT uses RF energy only for its internal functions. Therefore, its RF emissions are very low and unlikely to cause interference in nearby electronic devices.
- HUBER 360[®] FIIT may be used in all establishments, including domestic sites and sites that are directly connected to the low voltage public power grid, which supplies domestic buildings.
- HUBER 360[®] FIIT should be installed on a floor made of wood, concrete, or ceramic tile. If the floor is made of a synthetic material, relative humidity should be at least 30%.
- \bullet HUBER 360 $^{\circ}$ FIIT should be connected to a commercial- or hospital-grade power source.
- If the user needs the HUBER 360[®] FIIT to continue working even during power cuts, it should be connected to an uninterruptible power supply or to a battery.
- Portable and mobile communication devices that use radio frequencies should not be used in close proximity of any part of the HUBER 360[®] FIIT, including its cables.
- There may be interference near to devices labeled with the following symbol 柳

\rightarrow general warranty conditions

You have recently acquired an appliance distributed by LPG Systems or an LPG Systems approved distributor. It is the purchaser/ user's responsibility to find out from the local authorities the conditions and professional qualifications required before using the appliance.

The purchase of this equipment implies the legal acceptance by the purchaser/professional user of these general warranty conditions. If the appliance was sold to you by an approved LPG Systems distributor, the purchaser/ user should refer to the supplier's warranty conditions. These may in no way increase the undertakings made by LPG in these present warranty conditions.

The warranty can only be implemented and is only valid if the warranty slip has been duly filled out and returned to LPG Systems within two weeks of delivery, irrespective of the country. Warranty slips that are only partially completed will be rejected.

The appliance is guaranteed against manufacturing flaws and defects in the raw materials.

The warranty extends for the shorter of the following two periods: two (2) year OR four thousand (4000) hours of use from the invoice date. During this period, LPG Systems undertakes to exchange or repair free of charge, as quickly as possible, any part that LPG Systems acknowledges as defective, however LPG Systems does not undertake to replace the entire appliance. Traveling and living expenses for our technicians and transportation costs of the appliance or parts to and from the aftersales service workshop are not covered by this warranty.

Replacements and repairs performed within this warranty, with or without immobilization of the equipment, shall not have the effect of extending the warranty period. Replaced parts become the property of LPG Systems or the approved distributor.

No compensation shall be paid for loss of use.

Subject to other conditions hereafter, warranty shall apply if the purchaser/professional user has allowed LPG to proceed to necessary repair works.

Warranty exclusion:

- Damage occurring during transportation. Transportation of this equipment and/or spare parts is at the recipient's own risk.
 Before taking delivery, it is the recipient's responsibility to verify the state of the goods and to make a claim against the transport company in the manner usual in the delivery country.
- Non-observance of the installation and operating instructions, failure to carry out maintenance and/or negligence in maintaining the appliance, connection to a faulty electricity supply or a non-earthed electricity supply or a power supply whose voltage is different to the one indicated on the appliance.
- If an appliance is sold before the end of the warranty period, the warranty is transferred to the purchaser for the remaining warranty period, on the condition that:
- I The original invoice is provided.
- II That the initial vendor is informed of the sale.
- Modification, mounting of accessories or dismantling of the equipment.

- Any operation and/or intervention not specified in the LPG Systems Operator's Manual and performed on the equipment by the purchaser/user and/or any party not approved by LPG Systems.
- Use of consumables, spare parts, inappropriate components or parts not supplied by LPG Systems.
- Normal wear and tear of any of the equipment's parts resulting from normal usage.
- Damages or default resulting of accidental events (falls, impacts...)
 Damages or default resulting of natural disasters (lightening, water damages...), fire, negligence or abuses.

7

WARRANTY

\rightarrow GENERAL WARRANTY CONDITIONS

Failure to comply with the general warranty conditions during the warranty period and after its expiry may constitute an exonerating cause of liability of LPG Systems in case of damage attributable to the delivered products.

The purchaser/user is responsible for the use of the device and will assume full responsibility for any damages, including damages caused to third parties, resulting from the failure to observe the instructions for use of the device and / or resulting from an improper use.

Under no circumstances LPG SYSTEMS will be held liable for any intangible or indirect damages, including any commercial or financial loss, loss of profit, loss of earnings, and damage to the brand image.

LPG SYSTEMS' liability, for all causes (with the exception of personal injury) is limited to the amount of the defective device's price.

\rightarrow WARRANTY ACTIVATION

YOU CAN ACTIVATE YOUR WARRANTY ONLINE AT:

http://warranty.lpgsystems.com

\rightarrow OR FILL IN THE FORM AND RETURN TO LPG SYSTEMS[®]

Name:
Address:
Email:
HUBER 360® FIIT serial number:
Date:

Through the form above, LPG SYSTEMS, which is the data controller, operates a processing of your personal data, based on the sales agreement of your HUBER 360[®] FIIT device, in order to manage the guarantees that apply to it. LPG SYSTEMS only stores your data for the time required by the operations for which it was collected and in compliance with the regulations in force.

In accordance with the applicable regulations, you have the right to access, rectify and delete your personal data, a right of limitation or opposition to the treatment for legitimate reasons, a right of define general and specific guidelines on how you intend to exercise your rights after your death, the right to portability of your data, and the right to lodge a complaint with a supervisory authority.

You can exercise your rights by email at dataprivacy@lpgsystems.com or by mail at: LPG SYSTEMS, Customer Service Department, Technoparc de la Plaine, 30 rue Dr Abel, CS 90035 - 26000 VALENCE (FRANCE), with a photocopy of an identity document bearing your signature.

7

WARRANTY



HEADQUARTERS LPG SYSTEMS S.A.S.

TECHNOPARC DE LA PLAINE 30, RUE DU DR. ABEL - CS 900 35 - 26902 VALENCE CEDEX 09 - FRANCE TEL.: +33 (0)4 75 78 69 00 - FAX: +33 (0)4 75 42 80 85

INTERNATIONAL/MARKETING

2753, ROUTE DES DOLINES - BP 243 06905 SOPHIA-ANTIPOLIS CEDEX - FRANCE TEL.: +33 (0)4 92 38 39 00 - FAX: +33 (0)4 92 96 09 65

